

# TruPort Support Head

Think laterally




## Ideal conditions for anesthesia workstations with the horizontal TruPort Support Head

The reliable TruPort modularity combined with space saving horizontal orientation makes TruPort Horizontal the ideal solution for the anesthesia workstation. Positioned directly above the anesthesia machine, TruPort Horizontal provides more freedom of movement in the OR and enables efficient workflow around the anesthesia workstation. Through modular architecture, supply modules for gases, power and data can be simply repositioned or even extended at any time according to your requirements – even post installation.

## More flexibility for anesthesia

- 2 lengths available (500 mm / 19.7" and 750 mm / 29.5")
- To be combined with all TruPort arm systems
- Up to 4 side rails laterally mounted



Not all products/ options are available in all countries. For further information about our products or services, please contact your local Trumpf Medical representative or visit our webpage:

[www.trumpfmedical.com](http://www.trumpfmedical.com)

- Accessories can be attached with just one click
- Maximum pay load of the head is 80 kg / 176 lbs
- Suitable TruPort accessories available – matching your requirements, provided wirlessly and electrically

## More space in the OR

Compact architecture of the horizontal TruPort support head. Positioned directly below the ceiling allows maximum head-room and freedom of movement.

## Flexible solutions for the future

React quickly and efficiently to new requirements in the OR. The simple repositioning and possibility to retrofit supply modules for gases, power and data, TruPort Horizontal provides the flexible solution for the future.

## User-friendly control

User-friendly control handles and an optional remote control allow ergonomic control for the user. Compatible with all arm systems, TruPort Horizontal can be easily controlled in all situation.

Trumpf Medical reserves the right to make changes without notice in design, specifications and models. The only warranty Trumpf Medical makes is the express written warranty extended on the sale or rental of its products. ©2015 TRUMPF Medizin Systeme GmbH + Co. KG, ALL RIGHTS RESERVED.

TRUMPF Medizin Systeme GmbH + Co. KG  
Benzstraße 26 · 82178 Puchheim · Germany · Phone +49 89 80907-0 · E-Mail [info@trumpfmedical.com](mailto:info@trumpfmedical.com)

Flyer ID: 2017174, 30/09/2015