


TruSystem™ 7000dV OR Table

Redefining Patient Positioning in *da Vinci*® Surgery



The Trumpf Medical™ TruSystem™ 7000dV Surgical Table is designed to work with Integrated Table Motion for the *da Vinci*® Xi™ Surgical System from Intuitive Surgical.

Through an exclusive partnership, Intuitive Surgical and Trumpf Medical take operating room integration to a new level.

Integrated Table Motion

The *da Vinci*® Xi™ Surgical System featuring Integrated Table Motion connects to Trumpf Medical’s TruSystem™ 7000dV OR Table so that a patient can be dynamically positioned while the surgeon operates. This allows the surgical team to efficiently manage the patient and utilize gravity for optimal exposure and access to the target anatomy.

TruSystem™ 7000dV OR Table Package

The TruSystem™ 7000dV OR Table is the central piece of a unique package which provides components and accessories for multi-specialty *da Vinci*® Surgery:

- Double Joint Head Section
- Two Parts Leg Section, Spreadable
- dV Cable Remote Control
- dV Connection Cable
- Universal Arm Supports
- Patient Restraint Strap

With any TruSystem™ 7000dV OR Table two-years of comprehensive service coverage and maintenance are included.



Patient Weight	1,000 lbs /450 kg	Exceptional weight capacity in any articulation
Low Height	23“(±0.4”) /59 cm (±1 cm)	Low platform for docking access
Trendelenburg Tilt	± 45° ± 30°	Superior flexibility for reaching target anatomy

Additional Advanced Positioning Accessories

Trumpf Medical offers an additional accessory package of specialty components that optimize patient positioning capabilities for advanced and complex robotic-assisted surgeries.



Not all products / options are available in all countries. For further information about our products or services, please contact your local Trumpf Medical representative or visit our webpage:

www.trumpfmedical.com

Trumpf Medical reserves the right to make changes without notice in design, specifications and models. The only warranty Trumpf Medical makes is the express written warranty extended on the sale or rental of its products.

©2016 Trumpf Medizin Systeme GmbH + Co. KG, ALL RIGHTS RESERVED.

TRUMPF Medizin Systeme GmbH + Co. KG
Benzstraße 26 · 82178 Puchheim · Germany · Phone +49 89 80907-0 · E-Mail info@trumpfmedical.com

Flyer ID: 2067173, 13.10.2016

*TruSystem® s a registered trademark of TRUMPF GmbH + Co. KG

*da Vinci® is a registered trademark of Intuitive Surgical, Inc.

*Xi™ is a trademark of Intuitive Surgical, Inc.

Enhancing outcomes for patients and their caregivers:

Hill-Rom