

Compella™ Bariatric bed

Advancing Bariatric Care



Features and benefits

Maintains bariatric **patient dignity** by providing a care setting in a bed designed to resemble a standard hospital bed

Reduce caregiver strain and workload through **one-button powered width expansion** and **FlexAfoot™** allowing for quick configuration of the bed

Drive efficiency and ease of use with **SlideGuard™** designed to minimise the occurrence of the patient sliding down in bed

Ensure safe patient care with the optional **IntelliDrive™ XL Power Transport**, which offers caregivers the ability to move the bed and patient with ease

Optimise patient care with **head-of-bed alarm** and optional **3-mode bed exit alarm** that alerts caregivers of any safety concerns

The sturdy siderails, **low bed position combined with seat deflate** and ergonomic controls provide patients the confidence to safely enter and exit the bed with greater independence and dignity

The **in-bed scale** enables caregiver to monitor the patients' weight from the comfort of the bed

Enhancing outcomes for patients and their caregivers:

Hill-Rom

Compella Bariatric bed



One-button powered width expansion



Optional: IntelliDrive XL



Point-of-Care™ caregiver controls



Optional: Patient pendant control

Technical Specifications

Safe Working Load	500 kg
Patient weight	113-454 kg
Height adjustment	47-68 cm
Transport height	43 cm
Protection class	IPX4
Bed angles	
Head section	0-48°
Knee section	0-16°
Reverse Trendelenburg/ Trendelenburg	+9°/-9°
Sleep deck width	102-127 cm
Overall width	
minimum (width extenders retracted)	109 cm
maximum (width extenders extended)	135 cm
Sleep deck length	203-224 cm
Overall length	
minimum (foot section retracted)	231 cm
maximum (foot section extended)	256 cm

Surfaces

- non-powered foam
- powered air with Turn Assist
- powered air with Turn Assist and CLRT (Continuous Lateral Rotation Therapy)

Castors	4 x 150 mm double-band
---------	------------------------

Controls

standard	static patient and caregiver
optional	patient pendant

Additional options

- 3-mode bed exit alarm
- IntelliDrive XL

Hill-Rom is a leading global medical technology company with more than 7,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

CE These medical devices bear the CE mark and comply with all applicable regulations. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Medical beds/mattresses manufacturer:

Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom Industries SA, 188 Rue du Caducée, Parc EuroMédecine, 34195 Montpellier Cedex 5 - France

Völker GmbH, Wullener Feld 79, 58454 Witten - Germany

Hill-Rom Holdings, Inc., 1069 State Route 46 East, Batesville, IN 47006 - USA

Lifts manufacturer: Liko AB, Nedre vägen 100, 975 92 Luleå - Sweden

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN126305-01, 05 March 2015

Enhancing outcomes for patients and their caregivers:

Hill-Rom