

Transport chair - Anatome

Makes transfer easier, enhances comfort and safety



Features and benefits

- Adjusts to individual patients' shape, meaning improved comfort during transfers
- Intuitive design enhances patient comfort and is easy to operate for the caregiver
- Removable back and seat cushions allows for easier cleaning of the chair
- Height-adjustable armrests and retractable headrest enhance comfort and safety for patients
- Retractable footrest allows for a more comfortable seated position
- Multiple IV pole sockets and accessory hooks offer flexibility when placing accessories

Transport chair - Anatome



Head support



Comfort seating pad



Overchair table



IV pole

Technical Specifications

model SM647B

Height	1,500 mm
Seat height (hydraulic)	495 - 715 mm
Width (with armrest)	680 mm
Depth (with footrest)	1,090 mm
Weight	58 kg
Maximum load	135 kg

Material

Frame	Epoxy coated welded metal
Armrests	Polyurethane
Fabric	Coated polycotton blend

Castors

standard	4 x 125 mm single-band with central brake
optional	4 x 125 mm double-band with central brake

Fabrics

see coated fabric colours sheet

Accessories

- head support
- comfort seating pad
- overchair table
- telescopic IV pole (4 hooks)

Hill-Rom is a leading global medical technology company with more than 10,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

The Anatome chair for adult patients is designed for general care units, surgery, day surgery and medicalised long-term care. It is used to transport the patient or transfer a patient from the chair to the bed and vice versa.

Class I

Manufacturer: Hill-Rom S.A.S. – BP 14 Zi du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN150323-03, 25. August 2017

Enhancing outcomes for patients and their caregivers:

Hill-Rom